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CLOSED-EAR SENSATION

What is ear wax?

Ear wax is produced by the ceruminous and sebaceous glands in the ear canal. Its precise function is to defend and protect the ear canal, constantly lubricating and moistening it to preserve its delicate physiological balance.

The presence of wax in the ear canal helps to prevent the entry of external pollutants and dust.

Normally the wax is propelled toward the outside of the ear by imperceptible skin movements.

Sometimes wax can become a problem. Why is this?

When the sebaceous glands produce too much wax, or when there is insufficient lubrication of the ear canal, the natural cleaning mechanism becomes blocked, making it difficult or impossible for the wax to be expelled from the ear. Following immersion in water, the water makes the wax swell up, blocking the ear canal. This often leads to the formation of a wax "plug", causing extremely irritating symptoms such as loss of hearing, whistling and buzzing in the ears, pain and even dizziness. In these cases, the use of inappropriate devices, such as cotton buds, to remove the wax should be avoided since this actually risks making the situation worse by compacting the wax deeper down.

Who is most affected by these problems?

• Children

It often happens that children's ceruminous glands are highly active, producing excess wax that leads to the formation of wax plugs.

• The elderly

With increasing age the glands become less active and produce dryer wax that is less able to move towards the exterior and thus it tends to accumulate in the ear canal.

• Hearing-aid wearers and earphone users

People who wear hearing aids or regular earphone users have much greater difficulty keeping the ear canal clean and free from wax and other impurities.

• Swimmers and divers

Water getting into the ear canal can make the wax swell up until it causes a blockage.

• Air travelers or people spending time in the mountains

These people are subject to repeated changes in ear pressure that can cause a number of annoying symptoms, particularly when wax plugs are present.

How to prevent excess ear wax

Regular use of Otosan® Spray is effective in both removing excess wax from the ears and preventing its formation.

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
EN Non-invasive cleansing

Otosan[®] Spray with organic plant extracts delivers delicate but effective action for cleaner, healthier ears because it eliminates excess wax from the ear canal and, if used regularly, it prevents wax formation.

Highly effective ingredients

Otosan[®] Spray with organic plant extracts is free from artificial colourings and preservatives in order to guarantee maximum levels of tolerance. Combines the effective cleansing action of saline solution with the antiseptic properties of tea tree oil. The moisturising action of aloe combines with the other properties of the 100% organically grown plant extracts: Lemon with outstanding purifying properties, the effective soothing action of Camomile and Marigold and the emollient power of Mallow.


How to use

With your head straight, bring the Otosan[®] Spray up to the ear so that the spray nozzle is in contact with the ear canal. Press the spray button a couple of times keeping the spray container upright.

Wait a few minutes to allow Otosan[®] Spray to do its work. Lean your head over to the side of the treated ear to let the liquid and impurities flow out. Clean and dry the ear gently with a tissue.

Repeat the treatment on the other ear. Rinse the spray nozzle with warm water after each use.

When to use

Otosan[®] Spray is a system that helps to maintain a good ear hygiene and to prevent earwax plugs. We recommend using it once a day.

Warnings

Use at room temperature (20°C-25°C). Do not expose the bottle to sources of heat. Do not use if there is damage or injury to the ear drum or canal, any ear related medical conditions, after recent surgery or in the case of allergies to any of the ingredients. Keep out of the reach of children. Do not use on children under 30 months of age. For external use only. Avoid contact with the eyes. Do not swallow. Do not use on broken skin. Do not use for any other purpose and in different ways than those specified by the manufacturer in these directions. Not to be used after the expiry date.

Side effects

To date, there are no known side effects.

OTOSAN SPRAY, OTOSAN CONES AND OTOSAN DROPS
a complete range of products
for cleaner, healthier ears